

Innovation



For who you are & what you need...

eXcursion™

Introducing the **lightest side-fold rollator around!**

At just 15.4 lb the *Airgo® eXcursion* weighs 20% less than the competition while still providing peace of mind with industry leading stability.



Easy to fold



Ultra-light frame



Ergo-form handgrips



Collapsible basket

www.airgomobility.com

Innovation

Airgo[®] eXcursion[™]

For who you are & what you need...



Padded seat and backrest for added comfort

Locks in open position for added safety

Reflective trim for added visibility in the dark

Soft basket doesn't need to be removed when folding

Ergo-form handgrips fit comfortably and place your hand at a natural angle for better brake access

Easy to activate hand brakes

Easy to fold by just a pull of the strap

Lightweight frame makes lifting and storage simple

8" Comfort-Grip non-marring wheels are ideal for indoor & outdoor use

Lime



Colour	Item No.	Item UPC
Pearl Black	700-920	775757009209
Pacific Blue	700-921	775757009216
Cranberry	700-922	775757009223
Lime	700-923	775757009230



5' - 6' 2"	30.5" - 37"	20"	17.5"	25"	28.5"
153 - 188cm	77 - 94cm	51cm	44cm	63cm	72cm