

REDGUM BRAND
Daily Living and
Lifestyle Products



FOLDING REACHER

RGMF32

ASK FOR THE BRAND

redgumbrand.com.au

FOLDING REACHER

FEATURES

- Excellent aid to extend user's reach
- Rotatable head 0°-90°
- Folds for storage
- Minimal hand grip required to operate jaws
- Comfortable ergonomic soft rubber grip on handle
- Large multi finger trigger
- Rubber grip on jaws
- Magnet on top of jaw
- Lightweight aluminium arm with plastic parts

SPECIFICATIONS

Overall Dimensions:

L: 810mm x W: 40mm x D: 130mm

Folded:


L: 470mm x W: 80mm x D: 130mm

Jaw Width: 100mm

Product weight: 240g

Max. load weight: 5kg

Please contact **REDGUM Brand** if you require further information on this product



FOLDS FOR STORAGE

ROTATING HEAD

**SOFT RUBBER GRIP
LARGE TRIGGER**

IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use. The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

P: 08 9248 4180

E: orders@for-de.com.au

1 Business way Malaga Western Australia 6090

Part of the **For-de Group Pty Ltd**